


Deep learning-based analysis pipeline
for studies of biological contaminants

By: Avi Bajaj



Scope

- Bacterial/viral adhesion and resilience to surfaces is a global healthcare problem that was exacerbated by the (ongoing) COVID-19 pandemic.
- My group is exploring the possibility of surface topology modification as a physical mitigation strategy.



Scope – cont.

- We are currently evaluating polymer systems and topologies based on their adhesion resistance to biological contaminants.
- We are integrating a suite of [experimental \(microscopy-based\)](#) and [computational \(simulation-based\)](#) techniques.
- We need an efficient way to capture the trends from microscopy data without compromising on accuracy!



About me



- 2nd year PhD student
- Work regularly with Python and C/C++ to develop pre/postprocessing routines for various kinds of data
- Avid tennis player and lifter of weights

Key Learnings:

- Deep learning integration is a big step forward for scientific imaging, where image quality can be questionable
- Software can be GUI-reliant without being completely GUI-dependent
- It's up to the user to read the documentation and find a creative solution that suits their need(s)


Tools and Techniques

Software/packages:


- Fiji (ImageJ) v1.54f: <https://downloads.imagej.net/fiji/archive/20230710-2317/>
 - StarDist: <https://imagej.net/plugins/stardist>
- CSBDeep: <https://pypi.org/project/csbdeep/>
- Imageio: <https://pypi.org/project/imageio/>
- Tifffile: <https://pypi.org/project/tifffile/>
- NumPy: <https://pypi.org/project/numpy/>

Methods:

- Deep learning (cell/nuclei detection)
- Fluorescence microscopy



Source sample
images (*E. coli*)



Apply preprocessing
commands by running ImageJ
headlessly, i.e.,

- subtract background
- normalize brightness



Convert
postprocessed
images to 8-bit
(grayscale)



Convert
preprocessed
images to 8-bit
(grayscale)




Run StarDist2D with pretrained
model/parameters, i.e.,


- 2D versatile (fluorescent nuclei)



Export number of
regions of interest
(ROIs) per image
to an external file




Results – cont.



Results – cont.



Results – cont.



Results – cont.



Not so lucky...

We can train a model with our own images or refine the segmentation parameters.

We can also try changing the image brightness during preprocessing.

Performance



Start-to-finish
processing time
of **~7.5 s/image**

This is **~97.5%**
decrease in
elapsed time
with the manual
workflow that
required us to
launch ImageJ
(5+ min/image)

Outcomes

- Developed an efficient cell counting pipeline that can be run entirely from cmd
- Large image artifacts are screened out
- Pre/postprocessed images are saved for qualitative comparisons and sanity checks
- Counts still vary with brightness settings – no “one size fits all”
- Smaller image artifacts are still being counted, so we may need to play with the segmentation algorithm parameters
- 8-bit conversion can be clunky

